


How to make cats wear Balto Carpal Brace.



1, Stick 3M double sided tape onto the inside of the brace.



2, Cut the bandage into appropriate size, then press the non-sticky side of bandage onto the 3M tape; Or, you can do it the other way around by sticking the bandage onto the cat's limb first .



3, Put the brace on (refer to the Balto instructions on how to do that), and wrap Petflex.

*I also stuck 2 bits of Velcro on the front side of her arm where the braces Velcro wraps. The bits can help to hold the braces Velcro in place to add more stability.

*Since the brace is 'glued' to her limb, thus there's no need for it to get wrapped too tight.